

# TERF2 Antibody, FITC conjugated



PACO46532

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TERF2

**Uniprot**

Q15554

**Synonyms:**

Telomeric repeat-binding factor 2 (TTAGGG repeat-binding factor 2) (Telomeric DNA-binding protein), TERF2, TRBF2 TRF2

**Immunogen:**

Recombinant Human Telomeric repeat-binding factor 2 protein (78-238AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A