

TERF2 Antibody, Biotin conjugated



PACO46533

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TERF2

Uniprot

Q15554

Synonyms:

Telomeric repeat-binding factor 2 (TTAGGG repeat-binding factor 2) (Telomeric DNA-binding protein), TERF2, TRBF2 TRF2

Immunogen:

Recombinant Human Telomeric repeat-binding factor 2 protein (78-238AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A