

KYNU Antibody, Biotin conjugated



PACO46549

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the cleavage of L-kynurenine (L-Kyn) and L-3-hydroxykynurenine (L-3OHKyn) into anthranilic acid, (AA) and 3-hydroxyanthranilic acid, (3-OHAA), respectively. Has a preference for the L-3-hydroxy form. Also has cysteine-conjugate-beta-lyase activity.

Gene ID:

KYNU

Uniprot

Q16719

Synonyms:

Kynureninase (EC 3.7.1.3) (L-kynurenine hydrolase), KYNU

Immunogen:

Recombinant Human Kynureninase protein (99-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A