

BRD3 Antibody, HRP conjugated



PACO46559

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds hyperacetylated chromatin and plays a role in the regulation of transcription, probably by chromatin remodeling and interaction with transcription factors. Regulates transcription by promoting the binding of the transcription factor GATA1 to its targets. Regulates transcription of the CCND1 gene.

Gene ID:

BRD3

Uniprot

Q15059

Synonyms:

Bromodomain-containing protein 3 (RING3-like protein), BRD3, KIAA0043 RING3L

Immunogen:

Recombinant Human Bromodomain-containing protein 3 protein (356-506AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A