

# BRD3 Antibody, Biotin conjugated



PACO46561

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds hyperacetylated chromatin and plays a role in the regulation of transcription, probably by chromatin remodeling and interaction with transcription factors. Regulates transcription by promoting the binding of the transcription factor GATA1 to its targets. Regulates transcription of the CCND1 gene.

**Gene ID:**

BRD3

**Uniprot**

Q15059

**Synonyms:**

Bromodomain-containing protein 3 (RING3-like protein), BRD3, KIAA0043 RING3L

**Immunogen:**

Recombinant Human Bromodomain-containing protein 3 protein (356-506AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A