

PACO46562

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Activates the protein kinase activity of mTORC1, and thereby plays a role in the regulation of apoptosis. Stimulates the phosphorylation of S6K1 and EIF4EBP1 through activation of mTORC1 signaling. Has low intrinsic GTPase activity.

Gene ID:

RHEB

Uniprot

Q15382

Synonyms:

GTP-binding protein Rheb (Ras homolog enriched in brain), RHEB, RHEB2


Immunogen:

Recombinant Human GTP-binding protein Rheb protein (1-182AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO46562 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO46562 at dilution of 1:100.