

RIC1 Antibody, FITC conjugated



PACO46576

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The RIC1-RGP1 complex acts as a guanine nucleotide exchange factor (GEF), which activates RAB6A by exchanging bound GDP for free GTP and may thereby required for efficient fusion of endosome-derived vesicles with the Golgi compartment. The RIC1-RGP1 complex participates in the recycling of mannose-6-phosphate receptors. Required for phosphorylation and localization of GJA1.

Gene ID:

RIC1

Uniprot

Q4ADV7

Synonyms:

RAB6A-GEF complex partner protein 1 (Connexin-43-interacting protein of 150 kDa) (Protein RIC1 homolog), RIC1, CIP150 KIAA1432

Immunogen:

Recombinant Human RAB6A-GEF complex partner protein 1 protein (951-1233AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A