

# CRB2 Antibody, HRP conjugated



PACO46587

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May play a role in polarized cells morphogenesis.

**Gene ID:**

CRB2

**Uniprot**

Q5IJ48

**Synonyms:**

Protein crumbs homolog 2 (Crumbs-like protein 2), CRB2

**Immunogen:**

Recombinant Human Protein crumbs homolog 2 protein (253-483AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A