

ANKS6 Antibody



PACO46594

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Required for renal function.

Gene ID:

ANKS6

Uniprot

Q68DC2

Synonyms:

Ankyrin repeat and SAM domain-containing protein 6 (Ankyrin repeat domain-containing protein 14) (SamCystin) (Sterile alpha motif domain-containing protein 6) (SAM domain-containing protein 6), ANKS6, ANKRD14 PKDR1 SAMD6


Immunogen:

Recombinant Human Ankyrin repeat and SAM domain-containing protein 6 protein (1-132AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO46594 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO46594 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).