

Thrsp Antibody, FITC conjugated



PACO46612

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in the regulation of lipogenesis, especially in lactating mammary gland. Important for the biosynthesis of triglycerides with medium-length fatty acid, chains. May modulate lipogenesis by interacting with MID1IP1 and preventing its interaction with ACACA. May function as transcriptional coactivator. May modulate the transcription factor activity of THRB.

Gene ID:

Thrsp

Uniprot

Q62264

Synonyms:

Thyroid hormone-inducible hepatic protein (Spot 14 protein) (S14) (SPOT14), Thrsp, S14

Immunogen:

Recombinant Mouse Thyroid hormone-inducible hepatic protein (1-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A