

METTL21C Antibody, HRP conjugated



PACO46627

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Protein-lysine methyltransferase.

Gene ID:

METTL21C

Uniprot

Q5VZV1

Synonyms:

Protein-lysine methyltransferase METTL21C (EC 2.1.1) (Methyltransferase-like protein 21C), METTL21C, C13orf39

Immunogen:

Recombinant Human Protein-lysine methyltransferase METTL21C protein (6-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A