

NCR3LG1 Antibody, Biotin conjugated



PACO46641

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Triggers NCR3-dependent natural killer cell activation.

Gene ID:

NCR3LG1

Uniprot

Q68D85

Synonyms:

Natural cytotoxicity triggering receptor 3 ligand 1 (B7 homolog 6) (B7-H6), NCR3LG1, B7H6

Immunogen:

Recombinant Human Natural cytotoxicity triggering receptor 3 ligand 1 protein (25-262AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A