

# CLEC9A Antibody, Biotin conjugated



PACO46645

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Functions as an endocytic receptor on a small subset of myeloid cells specialized for the uptake and processing of material from dead cells. Recognizes filamentous form of actin in association with particular actin-binding domains of cytoskeletal proteins, including spectrin, exposed when cell membranes are damaged, and mediate the cross-presentation of dead-cell associated antigens in a Syk-dependent manner.

**Gene ID:**

CLEC9A

**Uniprot**

Q6UXN8

**Synonyms:**

C-type lectin domain family 9 member A, CLEC9A

**Immunogen:**

Recombinant Human C-type lectin domain family 9 member A protein (103-214AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A