

DARS2 Antibody, Biotin conjugated



PACO46665

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial matrix, mitochondrion, nucleus, aspartate-tRNA ligase activity, aspartate-tRNA(Asn) ligase activity, ATP binding, protein homodimerization activity, tRNA binding, mitochondrial asparaginyl-tRNA aminoacylation, tRNA aminoacylation.

Gene ID:

DARS2

Uniprot

Q6PI48

Synonyms:

Aspartate--tRNA ligase, mitochondrial (EC 6.1.1.12) (Aspartyl-tRNA synthetase) (AspRS), DARS2

Immunogen:

Recombinant Human Aspartate--tRNA ligase, mitochondrial protein (489-628AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A