

SLC30A10 Antibody, Biotin conjugated



PACO46669

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a pivotal role in manganese transport. Manganese is an essential cation for the function of several enzymes, including some crucially important for the metabolism of neurotransmitters and other neuronal metabolic pathways.

Gene ID:

SLC30A10

Uniprot

Q6XR72

Synonyms:

Zinc transporter 10 (ZnT-10) (Manganese transporter SLC30A10) (Solute carrier family 30 member 10), SLC30A10, ZNT10 ZNT8

Immunogen:

Recombinant Human Zinc transporter 10 protein (58-278AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A