

# WDFY4 Antibody, FITC conjugated



PACO46672

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

endomembrane system, extrinsic component of membrane, phospholipid binding.

**Gene ID:**

WDFY4

**Uniprot**

Q6ZS81

**Synonyms:**

WD repeat- and FYVE domain-containing protein 4, WDFY4, C10orf64 KIAA1607

**Immunogen:**

Recombinant Human WD repeat- and FYVE domain-containing protein 4 protein (1-322AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A