

WDFY4 Antibody, Biotin conjugated



PACO46673

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

endomembrane system, extrinsic component of membrane, phospholipid binding.

Gene ID:

WDFY4

Uniprot

Q6ZS81

Synonyms:

WD repeat- and FYVE domain-containing protein 4, WDFY4, C10orf64 KIAA1607

Immunogen:

Recombinant Human WD repeat- and FYVE domain-containing protein 4 protein (1-322AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A