

ZNF365 Antibody, FITC conjugated



PACO46696

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ZNF365

Uniprot

Q70YC5

Synonyms:

Protein ZNF365 (DISC1-binding zinc-finger protein) (Protein su48), ZNF365, DBZ KIAA0844

Immunogen:

Recombinant Human Protein ZNF365 protein (241-401AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A