

CARM1 Antibody, HRP conjugated



PACO46719

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CARM1

Uniprot

Q86X55

Synonyms:

Histone-arginine methyltransferase CARM1 (EC 2.1.1.319) (Coactivator-associated arginine methyltransferase 1) (Protein arginine N-methyltransferase 4), CARM1, PRMT4

Immunogen:

Recombinant Human Histone-arginine methyltransferase CARM1 protein (209-379AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A