

DAAM2 Antibody



PACO46734

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

extracellular exosome.

Gene ID:

DAAM2

Uniprot

Q86T65

Synonyms:

Disheveled-associated activator of morphogenesis 2, DAAM2, KIAA0381


Immunogen:

Recombinant Human Disheveled-associated activator of morphogenesis 2 protein (1-101AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO46734 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO46734 at dilution of 1:100.