

DAAM2 Antibody, Biotin conjugated



PACO46737

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular exosome.

Gene ID:

DAAM2

Uniprot

Q86T65

Synonyms:

Disheveled-associated activator of morphogenesis 2, DAAM2, KIAA0381

Immunogen:

Recombinant Human Disheveled-associated activator of morphogenesis 2 protein (1-101AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A