

# MILR1 Antibody, Biotin conjugated



PACO46745

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Immunoglobulin-like receptor which plays an inhibitory role in degranulation of mast cells. Negatively regulates IgE-mediated mast cell activation and suppresses the type I immediate hypersensitivity reaction.

**Gene ID:**

MILR1

**Uniprot**

Q7Z6M3

**Synonyms:**

Allergin-1 (Allergy inhibitory receptor 1) (Mast cell antigen 32) (MCA-32) (Mast cell immunoglobulin-like receptor 1), MILR1, C17orf60 MCA32

**Immunogen:**

Recombinant Human Allergin-1 protein (249-343AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A