

ZBTB46 Antibody, Biotin conjugated



PACO46749

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Function as a transcriptional repressor for PRDM1.

Gene ID:

ZBTB46

Uniprot

Q86UZ6

Synonyms:

Zinc finger and BTB domain-containing protein 46 (BTB-ZF protein expressed in effector lymphocytes) (BZEL) (BTB/POZ domain-containing protein 4) (Zinc finger protein 340), ZBTB46, BTBD4 ZNF340

Immunogen:

Recombinant Human Zinc finger and BTB domain-containing protein 46 protein (174-362AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A