

# SULF1 Antibody, FITC conjugated



PACO46784

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Exhibits arylsulfatase activity and highly specific endoglucosamine-6-sulfatase activity. It can remove sulfate from the C-6 position of glucosamine within specific subregions of intact heparin. Diminishes HSPG (heparan sulfate proteoglycans) sulfation, inhibits signaling by heparin-dependent growth factors, diminishes proliferation, and facilitates apoptosis in response to exogenous stimulation.

**Gene ID:**

SULF1

**Uniprot**

Q8IWU6

**Synonyms:**

Extracellular sulfatase Sulf-1 (hSulf-1) (EC 3.1.6), SULF1, KIAA1077

**Immunogen:**

Recombinant Human Extracellular sulfatase Sulf-1 protein (105-372AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A