

SULF1 Antibody, Biotin conjugated



PACO46785

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Exhibits arylsulfatase activity and highly specific endoglucosamine-6-sulfatase activity. It can remove sulfate from the C-6 position of glucosamine within specific subregions of intact heparin. Diminishes HSPG (heparan sulfate proteoglycans) sulfation, inhibits signaling by heparin-dependent growth factors, diminishes proliferation, and facilitates apoptosis in response to exogenous stimulation.

Gene ID:

SULF1

Uniprot

Q8IWU6

Synonyms:

Extracellular sulfatase Sulf-1 (hSulf-1) (EC 3.1.6), SULF1, KIAA1077

Immunogen:

Recombinant Human Extracellular sulfatase Sulf-1 protein (105-372AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A