

P3H2 Antibody, Biotin conjugated



PACO46793

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Shows prolyl 3-hydroxylase activity catalyzing the post-translational formation of 3-hydroxyproline in -Xaa-Pro-Gly-sequences in collagens, especially types II, IV and V.

Gene ID:

P3H2

Uniprot

Q8IVL5

Synonyms:

Prolyl 3-hydroxylase 2 (EC 1.14.11.7) (Leprecan-like protein 1) (Myxoid liposarcoma-associated protein 4), P3H2, LEPREL1 MLAT4

Immunogen:

Recombinant Human Prolyl 3-hydroxylase 2 protein (301-528AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A