

ABHD5 Antibody, Biotin conjugated



PACO46817

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Lysophosphatidic acid, acyltransferase which functions in phosphatidic acid, biosynthesis. May regulate the cellular storage of triacylglycerol through activation of the phospholipase PNPLA2. Involved in keratinocyte differentiation. Regulates lipid droplet fusion.

Gene ID:

ABHD5

Uniprot

Q8WTS1

Synonyms:

1-acylglycerol-3-phosphate O-acyltransferase ABHD5 (EC 2.3.1.51) (Abhydrolase domain-containing protein 5) (Lipid droplet-binding protein CGI-58), ABHD5, NCIE2

Immunogen:

Recombinant Human 1-acylglycerol-3-phosphate O-acyltransferase ABHD5 protein (2-349AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A