

SIRT1 Antibody



PACO46822

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

SIRT1

Uniprot

Q96EB6

Synonyms:

NAD-dependent protein deacetylase sirtuin-1 (hSIRT1) (EC 3.5.1.).


Immunogen:

Recombinant Human NAD-dependent protein deacetylase sirtuin-1 protein (2-747AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO46822 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO46822 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).