

SIRT1 Antibody, FITC conjugated



PACO46824

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

SIRT1

Uniprot

Q96EB6

Synonyms:

NAD-dependent protein deacetylase sirtuin-1 (hSIRT1) (EC 3.5.1.).

Immunogen:

Recombinant Human NAD-dependent protein deacetylase sirtuin-1 protein (2-747AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4.

Product Images

N/A

N/A