

Nlrp3 Antibody, Biotin conjugated



PACO46849

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Nlrp3

Uniprot

Q8R4B8

Synonyms:

NACHT, LRR and PYD domains-containing protein 3 (Cold autoinflammatory syndrome 1 protein homolog) (Cryopyrin) (Mast cell maturation-associated-inducible protein 1) (PYRIN-containing APAF1-like protein 1), Nlrp3, Cias1 Mmig1 Nalp3 Pypaf1

Immunogen:

Recombinant Mouse NACHT, LRR and PYD domains-containing protein 3 protein (1-153AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A