

LRR1 Antibody, Biotin conjugated



PACO46869

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May negatively regulate the 4-1BB-mediated signaling cascades which result in the activation of NK-kappaB and JNK1. Probable substrate recognition subunit of an ECS (Elongin BC-CUL2/5-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Gene ID:

LRR1

Uniprot

Q96L50

Synonyms:

Leucine-rich repeat protein 1 (4-1BB-mediated-signaling molecule) (4-1BBirr) (LRR-repeat protein 1) (LRR-1) (Peptidylprolyl isomerase-like 5), LRR1, PPIL5

Immunogen:

Recombinant Human Leucine-rich repeat protein 1 protein (18-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A