

# SMARCE1 Antibody, Biotin conjugated



PACO46893

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SMARCE1

**Uniprot**

Q969G3

**Synonyms:**

SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E member 1 (BRG1-associated factor 57) (BAF57), SMARCE1, BAF57

**Immunogen:**

Recombinant Human SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E member 1 protein (20-243AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A