

DNAH7 Antibody, Biotin conjugated



PACO46901

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Force generating protein of respiratory cilia. Produces force towards the minus ends of microtubules. Dynein has ATPase activity; the force-producing power stroke is thought to occur on release of ADP.

Gene ID:

DNAH7

Uniprot

Q8WXX0

Synonyms:

Dynein heavy chain 7, axonemal (Axonemal beta dynein heavy chain 7) (Ciliary dynein heavy chain 7) (Dynein heavy chain-like protein 2) (hDHC2), DNAH7, KIAA0944

Immunogen:

Recombinant Human Dynein heavy chain 7, axonemal protein (730-903AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A