

CDS1 Antibody, HRP conjugated



PACO46903

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Provides CDP-diacylglycerol, an important precursor for the synthesis of phosphatidylinositol (PtdIns), phosphatidylglycerol, and cardiolipin. Overexpression may amplify cellular signaling responses from cytokines. May also play an important role in the signal transduction mechanism of retina and neural cells.

Gene ID:

CDS1

Uniprot

Q92903

Synonyms:

Phosphatidate cytidyltransferase 1 (EC 2.7.7.41) (CDP-DAG synthase 1) (CDP-DG synthase 1) (CDP-diacylglycerol synthase 1) (CDS 1) (CDP-diglyceride pyrophosphorylase 1) (CDP-diglyceride synthase 1) (CTP: phosphatidate cytidyltransferase 1), CDS1, CDS

Immunogen:

Recombinant Human Phosphatidate cytidyltransferase 1 protein (6-85AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A