

# ATG16L2 Antibody, HRP conjugated



PACO46911

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in autophagy.

**Gene ID:**

ATG16L2

**Uniprot**

Q8NAA4

**Synonyms:**

Autophagy-related protein 16-2 (APG16-like 2) (WD repeat-containing protein 80), ATG16L2, WDR80

**Immunogen:**

Recombinant Human Autophagy-related protein 16-2 protein (262-348AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A