

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Plays a role in the regulation of lipogenesis, especially in lactating mammary gland. Important for the biosynthesis of triglycerides with medium-length fatty acid, chains. May modulate lipogenesis by interacting with MID1IP1 and preventing its interaction with ACACA. May function as transcriptional coactivator. May modulate the transcription factor activity of THRβ.

**Gene ID:**

THRSP

**Uniprot**

Q92748

**Synonyms:**

Thyroid hormone-inducible hepatic protein (Spot 14 protein) (S14) (SPOT14), THRSP

**Immunogen:**


Recombinant Human Thyroid hormone-inducible hepatic protein (41-144AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human lung cancer using PACO46930 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO46930 at dilution of 1:100.