

TET1 Antibody, Biotin conjugated



PACO46953

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

TET1

Uniprot

Q8NFU7

Synonyms:

Methylcytosine dioxygenase TET1 (EC 1.14.11. n2) (CXXC-type zinc finger protein 6) (Leukemia-associated protein with a CXXC domain) (Ten-eleven translocation 1 gene protein), TET1, CXXC6 KIAA1676 LCX

Immunogen:

Recombinant Human Methylcytosine dioxygenase TET1 protein (1700-1800AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A