

SLC16A10 Antibody, Biotin conjugated



PACO46969

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sodium-independent transporter that mediates the uptake of aromatic acid, . Can function as a net efflux pathway for aromatic amino acid, in the basolateral epithelial cells.

Gene ID:

SLC16A10

Uniprot

Q8TF71

Synonyms:

Monocarboxylate transporter 10 (MCT 10) (Aromatic amino acid, transporter 1) (Solute carrier family 16 member 10) (T-type amino acid, transporter 1), SLC16A10, MCT10 TAT1

Immunogen:

Recombinant Human Monocarboxylate transporter 10 protein (159-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A