

JARID2 Antibody, HRP conjugated



PACO46979

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

JARID2

Uniprot

Q92833

Synonyms:

Protein Jumonji (Jumonji/ARID domain-containing protein 2), JARID2, JMJ

Immunogen:

Recombinant Human Protein Jumonji protein (1-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A