

PRMT6 Antibody, FITC conjugated



PACO47008

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

PRMT6

Uniprot

Q96LA8

Synonyms:

Protein arginine N-methyltransferase 6 (EC 2.1.1.319) (Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 6) (Histone-arginine N-methyltransferase PRMT6), PRMT6, HRMT1L6

Immunogen:

Recombinant Human Protein arginine N-methyltransferase 6 protein (119-375AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A