

PEBP4 Antibody



PACO47010

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Seems to promote cellular resistance to TNF-induced apoptosis by inhibiting activation of the Raf-1/MEK/ERK pathway, JNK and phosphatidylethanolamine externalization.

Gene ID:

PEBP4

Uniprot

Q96S96

Synonyms:

Phosphatidylethanolamine-binding protein 4 (PEBP-4) (hPEBP4) (Protein cousin-of-RKIP 1), PEBP4, CORK1


Immunogen:

Recombinant Human Phosphatidylethanolamine-binding protein 4 protein (23-227AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO47010 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO47010 at dilution of 1:100.