

# SMAD5 Antibody, Biotin conjugated



PACO47021

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcriptional modulator activated by BMP (bone morphogenetic proteins) type 1 receptor kinase. SMAD5 is a receptor-regulated SMAD (R-SMAD).

**Gene ID:**

SMAD5

**Uniprot**

Q99717

**Synonyms:**

Mothers against decapentaplegic homolog 5 (MAD homolog 5) (Mothers against DPP homolog 5) (JV5-1) (SMAD family member 5) (SMAD 5) (Smad5) (hSmad5), SMAD5, MADH5

**Immunogen:**

Recombinant Human Mothers against decapentaplegic homolog 5 protein (1-209AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A