

WAC Antibody, Biotin conjugated



PACO47025

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a linker between gene transcription and histone H2B monoubiquitination at 'Lys-120' (H2BK120ub1). Interacts with the RNA polymerase II transcriptional machinery via its WW domain and with RNF20-RNF40 via its coiled coil region, thereby linking and regulating H2BK120ub1 and gene transcription. Regulates the cell-cycle checkpoint activation in response to DNA damage. Positive regulator of amino acid, starvation-induced autophagy. May negatively regulate the ubiquitin proteasome pathway.

Gene ID:

WAC

Uniprot

Q9BTA9

Synonyms:

WW domain-containing adapter protein with coiled-coil, WAC, KIAA1844

Immunogen:

Recombinant Human WW domain-containing adapter protein with coiled-coil protein (10-102AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A