

IFT46 Antibody, Biotin conjugated



PACO47077

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Forms part of a complex involved in intraflagellar transport (IFT), the bi-directional movement of particles required for the assembly, maintenance and functioning of primary cilia. May play a role in chondrocyte maturation and skeletogenesis.

Gene ID:

IFT46

Uniprot

Q9NQC8

Synonyms:

Intraflagellar transport protein 46 homolog, IFT46, C11orf2 C11orf60

Immunogen:

Recombinant Human Intraflagellar transport protein 46 homolog protein (1-304AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A