

ABCG8 Antibody, FITC conjugated



PACO47100

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transporter that appears to play an indispensable role in the selective transport of the dietary cholesterol in and out of the enterocytes and in the selective sterol excretion by the liver into bile.

Gene ID:

ABCG8

Uniprot

Q9H221

Synonyms:

ATP-binding cassette sub-family G member 8 (Sterolin-2), ABCG8

Immunogen:

Recombinant Human ATP-binding cassette sub-family G member 8 protein (176-413AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A