

CWC22 Antibody, Biotin conjugated



PACO47113

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for pre-mRNA splicing and for exon-junction complex (EJC) assembly. Hinders EIF4A3 from non-specifically binding RNA and escorts it to the splicing machinery to promote EJC assembly on mature mRNAs. Through its role in EJC assembly, required for nonsense-mediated mRNA decay.

Gene ID:

CWC22

Uniprot

Q9HCG8

Synonyms:

Pre-mRNA-splicing factor CWC22 homolog (Nucampholin homolog) (fSAPb), CWC22, KIAA1604 NCM

Immunogen:

Recombinant Human Pre-mRNA-splicing factor CWC22 homolog protein (360-623AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A