

EHD3 Antibody, Biotin conjugated



PACO47117

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

EHD3

Uniprot

Q9NZN3

Synonyms:

EH domain-containing protein 3 (PAST homolog 3), EHD3, EHD2 PAST3

Immunogen:

Recombinant Human EH domain-containing protein 3 protein (269-498AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A