

GMIP Antibody, HRP conjugated



PACO47119

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Stimulates, in vitro and in vivo, the GTPase activity of RhoA.

Gene ID:

GMIP

Uniprot

Q9P107

Synonyms:

GEM-interacting protein (GMIP), GMIP

Immunogen:

Recombinant Human GEM-interacting protein (719-858AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A