

HDHD2 Antibody



PACO47134

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

extracellular exosome, enzyme binding, phosphatase activity, dephosphorylation.

Gene ID:

HDHD2

Uniprot

Q9H0R4

Synonyms:

Haloacid, dehalogenase-like hydrolase domain-containing protein 2, HDHD2


Immunogen:

Recombinant Human Haloacid, dehalogenase-like hydrolase domain-containing protein 2 protein (97-249AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO47134 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO47134 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).