

FLRT3 Antibody, HRP conjugated



PACO47151

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

FLRT3

Uniprot

Q9NZU0

Synonyms:

Leucine-rich repeat transmembrane protein FLRT3 (Fibronectin-like domain-containing leucine-rich transmembrane protein 3), FLRT3, KIAA1469

Immunogen:

Recombinant Human Leucine-rich repeat transmembrane protein FLRT3 protein (209-264AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A